



# Ophthalmic Stand SN-500

*Serving Your Vision*



# Ophthalmic Stand SN-500

A versatile ophthalmic stand designed for reducing user's workload. A wide range of optical instruments can be set up according to user's preference.

## Outstanding features

### Overhead lamp

The overhead lamp knob on the control panel can dim up and down easily ensuring easy eye examination in the most suitable environment.

### View Tester arm

A wide range of vertical movement and outstanding counterbalancing ensure smooth and easy operation with a view tester mounted.

### Refractometer arm

A refractometer, etc. on an optional table top can move 220mm vertically.

### Lowest arm

A slit lamp on an optional table top can move up and down easily by the up/down button under the table top.

\* A power-supply voltage from each arm's power outlet is 50VA.

### Control box

The center and left rechargeable instruments wells can accommodate a rechargeable ophthalmoscope each (DC 4V)\*.

The right well can accommodate an ophthalmoscope connected to the binding post.

The output voltage of the binding post can be selected from 6 levels between AC 2.5V and AC 12V by the BD output voltage selection switch on the back.

\* Some ophthalmoscopes cannot be charged depending on these battery types or sizes.

\* Ophthalmoscopes with lithium ion batteries cannot be charged in the wells.

## Specifications

<b>Model name</b>	SN-500 Ophthalmic Stand
<b>Space required</b>	535 × 535 mm
<b>Weight</b>	137 kg




Operation panel



SN-500 combined with Ophthalmic Chair UN-15

\* VT-5 View tester, CP-40 Chart projector, UN-15 Ophthalmic chair, U08-04 Projector arm, and U12-05 Slit lamp table top are separately available.

 To ensure safe and correct handling please read the user manual before using.

- Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.
- Colours as they appear in brochure photographs may be slightly different to actual products due to lighting conditions when photographs were taken, or print colours.



**For The Americas, Asia-Pacific & Middle East**  
**Takagi Seiko Co., Ltd**

330-2 Iwafune, Nakano-shi, Nagano-ken, 383-8585, Japan  
TEL : +81(0)269-22-4511(Switchboard) URL : <https://www.takagi-j.com>



**For Europe & Africa**

**Takagi Ophthalmic Instruments Europe Ltd**

Citylabs 1.0, Nelson Street, Manchester, M13 9NQ, UK  
TEL : +44 (0)161 273 6330 URL : <https://www.takagieurope.com>



Your local distributor