

MS1-Compatible Smartphone Reference Diagram

Precautions for Use

- As a general guideline, smartphones compatible with the MS1 should be approximately 66-80 mm wide and approximately 138-160 mm long. Compatible models should also have the main camera positioned in the upper left as viewed from the rear.
- Models for which attachment has been confirmed by TAKAGI are shown below as “Verified Devices”.
- Other models confirmed to meet the dimensional requirements as of June 2026 are shown below as “Unverified Devices”.
- Please note that, even if the dimensional requirements are met, the smartphone may not be attachable depending on the shape, positions, and protrusion of the camera module.
- Even if the smartphone can be attached, it may be tilted depending on its shape. In this case, use spacers or similar items as necessary to adjust the smartphone so that it is level.
- This diagram is provided as a reference for determining whether a smartphone can be attached. Actual product and smartphone shapes may differ from those shown.

