

## AT-1 Applanation Tonometer

|                               |                                                               |      |
|-------------------------------|---------------------------------------------------------------|------|
| Model name                    | AT-1 Applanation Tonometer                                    |      |
| Attachable slit lamp models   | 700GL, 700GL NSW, 700GL LPH, 4ZL, 2ZL                         | 30GL |
| Mounting tonometer base       | TB700                                                         | TB30 |
| Measurement range             | 0-80 mmHg                                                     |      |
| Minimum scale value           | 2 mmHg                                                        |      |
| Measurement method            | Goldmann method (Goldmann Applanation Tonometry)              |      |
| Measurement deviation         | Nominal value $\pm 1.5\%$ or $\pm 0.49$ mN, whichever greater |      |
| Reverse span (Hysteresis)     | $\leq 0.29$ mN                                                |      |
| Diameter of applanation field | 3.06 mm                                                       |      |
| Depression area               | 7.354 mm <sup>2</sup>                                         |      |
| Dimensions                    | 72 (W) x 135 (D) x 251 (H) mm                                 |      |
| Weight                        | 0.44 kg                                                       |      |

\*Specifications as detailed in this sheet may be subject to change as improvements are made to products.

\*For the latest information, please visit our website (<https://www.takagieurope.com/>).