

AT-2 Applanation Tonometer

Model name	AT-2 Applanation Tonometer
Attachable slit lamp models	4ZL, 2ZL
Mounting tonometer base	TB2
Measurement range	0-80 mmHg
Minimum scale value	2 mmHg
Measurement method	Goldmann method (Goldmann Applanation Tonometry)
Measurement deviation	Nominal value $\pm 1.5\%$ or ± 0.49 mN, whichever greater
Reverse span (Hysteresis)	≤ 0.29 mN
Diameter of applanation field	3.06 mm
Depression area	7.354 mm ²
Dimensions	84.3 (W) x 207.6 (D) x 200.6 (H) mm
Weight	0.67 kg

*Specifications as detailed in this sheet may be subject to change as improvements are made to products.

*For the latest information, please visit our website (<https://www.takagieurope.com/>).