

CR-750/CR-750S Motorized Patient Chair

Model name	CR-750 Motorized Patient Chair
	CR-750S Motorized Patient Chair
Seat height	400-550 mm
Elevation stroke	150 mm
Load capacity	150 kg
Input power	AC 100-230V, 50-60 Hz
Power consumption	120 VA
Duty cycle	10% (2 minutes of continuous use / 18 minutes of no use)
Dimensions	430 (W) x 439 (D) x 715-865 (H) mm
Weight	26.5 kg (CR-750)
	25 kg (CR-750S)
Colour variation	The America, Asia-Pacific & Middle East: Pink, Blue and Green Europe & Africa: Blue/ROHS2 compliant *Pink and Green are not available in Europe & Africa

*Specifications as detailed in this sheet may be subject to change as improvements are made to products.

*For the latest information, please visit our website (<https://www.takagieurope.com/>).