

MS1 Mobile Slit Lamp Microscope

Model name	MS1 Mobile Slit Lamp
Slit illumination angle	30° fixed from left
Slit width, length	0.2mm, 14 mm
Illumination field diameter	14mm
Filter	Blue filter
Light-intensity control	Continuously variable, ON/OFF
Power supply	Two AAA batteries
Weight	230 g (excluding a smartphone and batteries)

*Specifications as detailed in this sheet may be subject to change as improvements are made to products.

*For the latest information, please visit our website (<https://www.takagieurope.com/>).