

OC-1A / OC-1B



Hydraulic doctor chair available with or without arm rests

Hard and thick cushion

Since the seat uses a hard and thick cushioning material, there is little sinking of the body, which helps with reducing tiredness even after long-term use.

Special seat shape


The seat is of unique triangle shape which reduces the pressure on the rear thighs.

Backrest to support the waist

The backrest, which has been devised in the width and curvature that holds the upper part of the waist, not only receives the back, but also has the effect of maintaining a posture that can reduce the burden on the operator's waist.


Simple elevating section

The foot-operated lift unit inherits the simple and highly durable mechanism from the previous model, and the foot pedal allows to perform easy and reliable lifting/lowering/swing lock (seat rotation lock).



Dimensions

Dimensions (mm)



*The OC-1B does not include the armrest.

Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.

Colours of product(s) in brochure may slightly differ from actual product(s) due to lighting situation when photographed or print conditions.

 To ensure safe and correct handling please read the user manual before using



Armrest with cushion

The OC-1A armrest with cushion is slightly inclined inwards so that the arms are placed naturally.

Caster with big diameter and wide width

Casters with big diameter and wide width reduce the resistance and enable the operator to move easily.

Stable base section

Five arm shafts and a low centre of gravity hydraulic mechanism form the stable base section.

Specifications

Model name	OC-1A Hydraulic Medical Chair OC-1B Hydraulic Medical Chair
Slide stroke	Pedal, hydraulic type
Allowable load	187 mm
Power	380 (W) × 390 (D) mm
Power consumption	Up/down 170 mm
Weight	31 kg (OC-1A) 25 kg (OC-1B)

Takagi Ophthalmic Instruments Europe Ltd.

Citylabs 1.0
(Former Royal Eye Hospital Building)
Nelson Street
Manchester, M13 9NQ
United Kingdom
T: +44 (0)161 273 6330
W: <https://www.takagieurope.com>



Takagi Seiko Co., Ltd.



learn more



takagieurope.com

Your local distributor