

SN-500



Ophthalmic stand complete with arms and table

Overhead lamp

The overhead lamp knob on the control panel can dim up and down easily ensuring easy eye examination in the most suitable environment.

View Tester arm

A wide range of vertical movement and outstanding counterbalancing ensure smooth and easy operation with a view tester mounted.

Refractometer arm

A refractometer, etc. on an optional table top can move 220mm vertically.

Lowest arm

A slit lamp on an optional table top can move up and down easily by the up/down button under the table top.

* A power-supply voltage from each arm's power outlet is 50VA.

Control box

The centre and left rechargeable instruments wells can accommodate a rechargeable ophthalmoscope each (DC 4V)*. The right well can accommodate an ophthalmoscope connected to the binding post.

The output voltage of the binding post can be selected from 6 levels between AC 2.5V and AC 12V by the BD output voltage selection switch on the back.

* Some ophthalmoscopes cannot be charged depending on the battery types or sizes.

* Ophthalmoscopes with lithium ion batteries cannot be charged in the wells.

Specifications

Model name	SN-500 Ophthalmic Stand
Space required	535 × 535 mm
Weight	137 kg

Description and appearance as detailed in this brochure may be subject to change as improvements are made to products.

Colours of product(s) in brochure may slightly differ from actual product(s) due to lighting situation when photographed or print conditions.

 To ensure safe and correct handling please read the user manual before using



Operation panel



SN-500 combined with UN-15 Ophthalmic Chair

* VT-5 View tester, CP-40 Chart projector, UN-15 Ophthalmic chair, U08-04 Projector arm, and U12-05 Slit lamp table top are separately available.

Takagi Ophthalmic Instruments Europe Ltd.

Citylabs 1.0
(Former Royal Eye Hospital Building)
Nelson Street
Manchester, M13 9NQ
United Kingdom
T: +44 (0)161 273 6330
W: <https://www.takagieurope.com>



Takagi Seiko Co., Ltd.



learn more



[takagieurope.com](https://www.takagieurope.com)

Your local distributor